


### 1. Introduction

We present a study investigating the effect of deletion of the epigenetic regulation enzyme O-6-methyllysine transferase (OGT) on genome-wide 5mC and 5hmC levels in mouse embryonic stem cells (mESC) using a novel 6-base sequencing technology.

The O-6-methyllysine transferase OGT interacts robustly with all three mammalian Ten-Eleven Translocation (TET) methylcytosine dioxygenases. Using duet evoC 6-base sequencing (enabling individual discrimination of A, C, T, G, 5mC, and 5hmC), we show that deletion of the Ogt gene in mESC results in a widespread increase in the TET product 5hmC in both euchromatic and heterochromatic compartments, with concomitant reduction of the TET substrate 5mC at the same genomic regions. mESC treated with an OGT inhibitor also displayed increased 5hmC demonstrating OGT enzymatic activity is needed to suppress TET activity. This indicates that OGT restrains TET activity and limits outward DNA demethylation in a manner that requires the TET-OGT interaction and the catalytic activity of OGT. DNA hypomethylation in OGT-deficient cells was accompanied by deuppression of transposable elements (TEs) predominantly located in heterochromatin.


We suggest that OGT protects the genome against TET-mediated DNA demethylation and loss of heterochromatin integrity, preventing the aberrant increase in TE expression noted in cancer, autoimmune/inflammatory disease, cellular senescence and aging.

### 2. duet evoC 6-base sequencing [A,C,T,G,5mC,5hmC]




**Figure 1.** (a) TET-mediated demethylation pathway. (b) Duet multibioscience solution evoC – a 6-base sequencing technology that reads all four canonical bases plus 5mC and 5hmC via strand copy, 5mC copy and 5mC + 5hmC protection enzymatic steps. (c) The duet multibioscience solution evoC works as an end-to-end solution comprising reagents & bioinformatics pipeline.

### 3. OGT deletion decreases 5mC and increases 5hmC across the genome



**Figure 2.** (a) Left: Generation of Ogt KO mESC. Right: Western blots showing loss of OGT protein (top) and O-6-methyllysine modification (bottom) in whole-cell lysates. (b) 6-base sequencing shows increased 5hmC (top) and decreased 5mC (bottom) in Ogt KO compared to Ctrl Ogt mESC at whole genome resolution (left) or chromatin compartments (right). (c) Genomic browser view of 5hmC (top) and 5mC (bottom) in overlaid tracks of Ogt KO (orange) and Ctrl Ogt mESC. 10 kb windows averaged. Euchromatin (green), heterochromatin (blue). (d) Violin diagrams showing reciprocal gain of 5hmC (red) and loss of 5mC (blue) in Ogt KO mESC vs Ctrl Ogt mESC for individual CpGs (left) or 10kb windows (right). (e) Genomic browser view showing reciprocal gain of 5hmC (red) and loss of 5mC (blue) in Ogt KO vs Ctrl Ogt mESC. Plotted log2 ratios (Ogt KO/ Ctrl Ogt mESC) for 5hmC (top) and 5mC (bottom). (f) Increase of 5hmC after OGT inhibition with OSMI-4. 5hmC measured by flow cytometry after 4 days of OSMI-4 treatment. Bar graphs represent the mean  $\pm$  s.d. from three independent experiments. (g) Violin plots showing increased 5hmC levels (left) but unaltered 5mC levels (right - likely due to earlier growth arrest) in OSMI-4-treated versus control mESC cells.


### 3. Decreased 5mC and increased 5hmC result in increased TE expression



**Figure 3.** (a) Volcano plot of RNA-seq data showing differential TE expression in Ogt KO vs Ctrl Ogt mESC. Red dots: transcripts exhibiting  $\geq$  twofold expression change. (b) Violin plots of Ogt KO vs Ctrl Ogt mESC cells analysed for differential expression of uniquely mapped TEs. No. of upregulated (red) or downregulated (grey) TEs. (c) Violin plots of TE subfamilies in Ogt KO vs Ctrl Ogt mESC cells. (d) Violin plots of TE subfamilies in Ogt KO vs Ctrl Ogt mESC cells. Difference of median 5hmC levels between conditions (Ogt KO = Ctrl) displayed on top.

### 4. Conclusion

We report here that OGT globally restrains TET enzymatic activity and maintains DNA methylation genome-wide, in a manner that depends on OGT catalytic activity and the TET-OGT interaction (illustrated below). TET activity was unleashed upon Ogt gene deletion on OGT inhibition, resulting in ongoing demethylation defined by parallel increases in 5hmC and decreases in 5mC at overlapping CpGs, observed using duet multibioscience solution evoC 6-base sequencing. The discovery (in press) reveals a novel mechanism for maintaining genomic stability, with important implications for both development and disease.



### 5. References

- Fullgrabe J, et al. Simultaneous sequencing of genetic and epigenetic bases in DNA. *Nat Biotechnol*. 2023 Oct;41(10):1454-1464.
- Sepulveda H, et al. OGT prevents DNA demethylation and suppresses the expression of transposable elements in heterochromatin by restraining TET activity genome-wide. *Nat. Struct. Mol. Biol.* (2025) - in press

