

Refining liquid biopsy with the 6-base genome: Generating more information from cell free DNA


Fabio Puddu, Tom Charlesworth, Robert Crawford, Nick Harding, Riccha Sethi, Jamie Scotcher, Annelie Johansson, Ermira Lleshi, Aurelie Modat, Michael S Wilson, Páidí Creed

1. Introduction

Liquid biopsy for profiling of cell free DNA (cfDNA) in blood holds huge promise to transform how we experience and manage cancer by early detection and identification of residual disease and subtype. However, a standard blood draw yields an average of only 10 ng of cfDNA, of which DNA derived from the tumour is a small minority.

Genetic and methylation data together have been shown to be more powerful for the detection of early cancer than either alone. Constrained to measuring four states of information, existing NGS-based technologies sacrifice genetic information for methylation calling.


duet multimodal solution *evoC* is a new sequencing technology that simultaneously derives all four genetic bases without ambiguity in C or T calls alongside distinguishing 5-methylcytosine and 5-hydroxymethylcytosine (6-base data) in a single read from a single DNA molecule [1]. The technology consists of pre-sequencing library prep and post-sequencing analysis pipeline, providing single-base resolution of genetics and epigenetics at high accuracy




State	Standard Sequencing Protocol	Protocol with C→T deamination
1	A	A
2	C/mC/hmC	mC/hmC
3	G	G
4	T	C/T

4. Accurate detection of genetic variants in a CRC patient cohort

cfDNA is thought to enter the bloodstream through apoptosis or necrosis, with cfDNA from healthy and cancer tissues released into the blood of cancer patients. To assess the ability of *duet evoC* to generate multimodal information from liquid biopsy samples, we obtained and sequenced cfDNA from 87 individuals, ranging from healthy volunteers to stage IV colorectal cancer (CRC) patients. (A) Somatic variants were called using Mutect2 on *duet evoC* in tumour-only mode and the presence of pathogenic or likely pathogenic variants associated with CRC is shown. An increase in the prevalence of these variants from stage I to stage IV patients mirrors an increase in the amounts of ctDNA (B), as estimated by ichorCNA [3].




6. Copy Number Variation




Copy number information extracted with CNVkit [6] from 10 cfDNA samples for each stage shows an increasing prevalence of chromosomal aberrations in later stages of CRC. Note that some chromosomes appear to be consistently rearranged in several samples (blue arrows)

2. duet multimodal solution *evoC*




5. Accurate epigenetic information for CRC detection

(A) DMRs identified from stage IV CRC tissue in the TCGA-COAD cohort [4] using Infinium Human Methylation 450K arrays were reproduced from 5mC levels in cfDNA obtained from stage IV CRC patients using *duet evoC* (12 CRC; 24 healthy volunteers). Cell-free DMRs where 5mC was greater in stage IV were defined as hypermethylated, and vice versa ($p < 0.05$; t-test). (B) Within these DMRs from the TCGA-COAD cohort, *duet evoC* captured differential (hydroxy)methylation in cfDNA from stage IV CRC patients. (C) Analysing cfDNA with *duet evoC* produces 6-base readouts that can be summarised across regions with ease using *modality*, part of the *duet* analysis suite. (D) ROC curves demonstrating that combining 5mC and 5hmC features in cfDNA improved separation of stage I CRC patients from healthy volunteers, when compared with using 5mC or 5hmC alone. Candidate features were defined as 5mC and 5hmC fractions of regions defined as stage IV tissue DMRs in TCGA-COAD, calculated from cfDNA from 24 stage I CRC patients and 25 healthy volunteers. Using *glmnet* [5], generalised linear models were trained on either 5mC or 5hmC features, or both (5mC + 5hmC), and evaluated using a leave-one-out cross-validation (LOOCV) approach. (E) Candidate features ranked by the number of times they were selected in each train/test split during LOOCV of the (5mC + 5hmC) model. (F) Analysing 5hmC improved the ability to distinguish between CRC stages in genomic regions with subtle 5mC differences. 5hmC (but not 5mC) fractions of regions overlapping KIF3 and CDH4 were selected as features in every LOOCV split. Averages were taken across the stage I CRC patients and healthy volunteers from (D), and cfDNA from 12 additional stage IV CRC patients.




3. Accurate genetic and epigenetic data

Complete accurate genome and methylome information from *duet evoC*. Genomic information is provided in full at higher or equivalent accuracy as other genetic and epigenetic methods (Figure 3A).




Simultaneously full methylome information is provided at high accuracy for both 5mC and 5hmC (Figure 3B). Figure (A) uses data generated on the Genome-in-a-bottle reference materials [2]




9. Conclusions

We have presented data illustrating the potential of *duet evoC* for liquid biopsy. With *duet evoC* it is possible to obtain multi-modal information, including SNPs, methylation, hydroxymethylation, fragmentomics, copy number variation and novel 2 dimensional biomarkers, from a single low-input sample of cfDNA. Further we have demonstrated improved ability to differentiate between healthy and stage IV CRC using combined methylation and hydroxymethylation features. *duet evoC*, which provides 6-base sequencing, is available to order now as a product from biomodal.

7. Fragmentomics information



8. Variant-associated methylation



10. References

1. Simultaneous sequencing of genetic and epigenetic bases in DNA, Füllgrabe and Gosal et al., Nature Biotechnology (2023) (*duet multimodal solution technology paper*)
2. Zook et al., Extensive sequencing of seven human genomes to characterize benchmark reference materials. *Sci Data* 3 (2016)
3. Adalsteinsson, V.A., Ha, G., Freeman, S.S. et al. Scalable whole-exome sequencing of cell-free DNA reveals high concordance with metastatic tumors. *Nat Commun* 8, 1324 (2017)
4. TCGA-COAD, <https://portal.gdc.cancer.gov/projects/TCGA-COAD>
5. Tay JK, Narasimhan B, Hastie T. "Elastic Net Regularization Paths for All Generalized Linear Models." *Journal of Statistical Software*, 106(1), 1–31 (2023)
6. Talevich et al., CNVkit: Genome-Wide Copy Number Detection and Visualization from Targeted DNA Sequencing. *PLOS Computational Biology* (2016)