

Genetic and Epigenetic Study of Formalin-damaged DNA


Robert Crawford*, Jamie Scotcher, Fabio Puddu, Dan Brudzewsky, Jane D Hayward, Andrada Tomoni, David Morley, Páidí Creed, Philippa Burns, and Joanna D Holbrook

Introduction


DNA comprises molecular information stored in genetic and epigenetic bases, both of which are vital to our understanding of biology. **duet multiomics solution +modC** is a technology which sequences at base resolution the complete genetic sequence integrated with modified cytosine (modC). **duet multiomics solution +modC** has previously been demonstrated using genomic and cell free DNA. Here we demonstrate its performance relative to formalin-compromised DNA standards [1] (fcDNA) as well as DNA extracted from formalin-fixed and paraffin embedded (FFPE) samples of two colorectal cancer (CRC) patients.

FFPE samples represent an important resource for studying genetic and epigenetic information from archived tissues. However, DNA damage induced by formalin fixation (e.g. deamination, fragmentation or nucleic acid cross-linking) can lead to decreased data quality from next generation sequencing (NGS) relative to 'gold standard' fresh-frozen samples. Damage may manifest in lower library yields and insert sizes as well as higher duplication and lower coverage rates[2]. When studying formalin-damaged DNA with **duet multiomics solution +modC** we found expected changes in yield and insert size, as well as C:G>T:A mutations. Background modC levels decreased ~5% in CpG contexts and increased ~0.5% in CHH/CHG contexts between untreated and severe formalin damaged (DIN ≤ 2.0) standards. Importantly, no substantial differences in variant allele frequency (VAF) were observed for a set of reference genes, even at severe formalin damage. Overall, **duet multiomics solution +modC** is compatible with formalin damaged samples such as FFPE.


1. Using duet multiomics solution +modC libraries with fcDNA




2. Formalin damaged DNA is compatible with +modC




3. Nucleotide bias is largely unaffected by formalin



4. Accurate CpG modC calls in formalin-damaged DNA



5. Accurate VAF estimates in formalin-damaged DNA




6. duet multiomics solution +modC libraries with FFPE samples

NA12878: undamaged gDNA
FFPE1: mild-moderate damage (DIN=5.0)
FFPE2: moderate-severe damage (DIN=3.2)

duet multiomics solution +modC applied to FFPE samples. To determine if *duet multiomics solution +modC* is compatible with formalin-fixed paraffin-embedded samples, we used as input DNA extracted from FFPE samples derived from two colorectal cancer (CRC) patients. DNA extraction was performed using the chemagic truXTAC DNA FFPE Kit. FFPE1: Female, CRC Grade 1, Stage IIIB, age 70-85, DIN = 5.0. FFPE2: Male, CRC Grade 1, Stage I, age 70-85, DIN = 3.2. Processing extracted DNA followed the same pre-sequencing workflow and data processing described in Panels 1 B & D. Here, 80ng of sheared extracted DNA from each FFPE sample were processed in triplicate. As with fcDNA lower yield and insert size was observed for formalin-damaged samples relative to undamaged NA12878 gDNA controls.

7. FFPE-samples are compatible with +modC



Conclusion

In conclusion we demonstrate the compatibility of **duet multiomics solution +modC** with formalin-damaged DNA. While damage resulted in lower library yields and insert sizes relative to undamaged controls, high accuracy genetic and epigenetic base resolution data is produced even at severe levels of formalin-induced damage. No effect on allele frequency detection was observed for the fcDNA QMRs standards. We note small changes in modC calling dependent on sequence context: a decrease for CpG, and an increase for CHH/CHG. We hypothesise this may be due to slightly different formalin-induced deamination pathways acting on unmodified cytosines (predominately found at CHH/CHG) and modified cytosines (much more prevalent in CpG contexts).

- Quantitative Multiplex Reference Standard (Horizon Discovery): <https://horizondiscovery.com/en/reference-standards/products/quantitative-multiplex-reference-standard-gdna>
- Hedgegaard, Jakob et al. "Next-generation sequencing of RNA and DNA isolated from paired fresh-frozen and formalin-fixed paraffin-embedded samples of human cancer and normal tissue." *PLoS one* vol. 9,5 e98187. 30 May. 2014.